



Transfusion Medicine Solution



When your patients
are relying on you,
rely on us...

QuidelOrtho transfusion solutions are designed to give laboratories the tools they need to confidently deliver efficient and reliable results, so that patients can get the transfusions they need.

With more than 80 years of trusted partnership with transfusion specialists and over 6,500 instruments, every second a patient is safely transfused due to our pretransfusion testing.



Ortho Vision® Swift Analyzer

The Ortho Vision Swift analyzer offers the most diverse testing capabilities available today, including serial dilutions, selected cell panels and open channels allowing labs to automate more tests and realize improved workflow efficiency.

When your lab is running at optimal efficiency, you can make reliable, accurate results available faster for patients, with 95% of results becoming available in less than 39 minutes.*

The standardization of the Ortho Vision® Swift platform solution across the lab and networks makes it the #1 choice for hospitals with more than 6,500 instruments installed worldwide.



Ortho® Workstation & Ortho Optix® Reader

The Ortho Optix Reader, when paired with the Ortho Workstation, provides a safe, cost-effective solution that helps labs continue to deliver consistent results with the same confidence you would get with an automated system.

The high-resolution camera captures colour, grayscale, and front and back images of each test without operator intervention. The solution provides objective data for automated results interpretation, eliminating the guess work, reducing risk of error, and delivering standardised processes.



Ortho BioVue® Technology

The unique column agglutination technology of the Ortho BioVue system works simply and effectively, delivering fast, accurate and standardized results.

The minimal centrifugation time of 5 minutes, offers one of the shortest turnaround times available and ensures that a blood group can be made available in 8.5 minutes.

The comprehensive cassette and reagent menu offers the ability to automate 99% of testing needs.

e-Connectivity® Technology

e-Connectivity technology, delivers increased uptime whereby transfusion specialists can expect to experience 24% less service calls, faster resolution through proactive analysis and the provision of a dashboard that offers key performance data related to your laboratory.

The e-Connectivity dashboard allows more control, providing a clear overview of key metrics, including platform turnaround time, workload, and usage patterns.

Laboratory operations are complex. Our service expertise and operational tools ensure that your solution works in sync so that you can provide the best patient care without disruption.

Ortho Connect®

Ortho Connect is an easy-to-use middleware that provides a centralised dashboard, where transfusion specialists can review and approve results from anywhere and ensure a seamless flow of information across your laboratory network.

This powerful software helps control data and metrics, optimize resources, and elevate patient care within the lab and across the network — all from a single connection.



Solution Benefits

- Shortest turnaround times on the market thanks to the minimal centrifugation time of 5 minutes ensures that a blood group can be made available in 8.5 minutes.
- 95% of results are available in less than 39 minutes.
- Comprehensive cassette and reagent menu offers the ability to automate 99% of testing needs.
- Experience 24% less service calls through proprietary e-Connectivity technology.



At QuidelOrtho, we transform diagnostic data into answers, understanding and action, illuminating the path forward for all. For more than 80 years, we've pursued the unknown with a passion and purpose to improve health. And we'll continue to transform the power of diagnostics into a healthier future for all.

quidelortho.com

© QuidelOrtho 2023

PR-14298